● 一次性使用医用口罩 小包装 中文版



● 一次性使用医用口罩 小包装 英文版



● 一次性使用医用口罩 外箱包装 中文版



Yu'an (Henan) Holding Co., Ltd

ADD. No. 168, industrial park, zhangsanzhai Town, Changyuan City, Xinxiang City, Henan Province

| Part |

● 一次性使用医用口罩 外箱包装 英文版



● 医用外科口罩 小包装 中文版



ADD. No. 168, industrial park, zhangsanzhai Town, Changyuan City, Xinxiang City, Henan Province

● 医用外科口罩 外箱包装 中文版



0222:宇安邦股_外科口罩纸箱1(63cm×45cm×500cm) 2+0专蓝(参考调色: 梁菡C96M60Y15; 浅菡C50M18)



0330:宇安萨版_JL差外科口罩網箱(63cm×45cm×50cm) 2+0专語(参考網色:環際C100M60:粉語C46M16)



Certificate

ICR Polska/CE/V/ M3100135

Name and address of Yu An (Henan) Holding Limited Company

certificate owner: No. 168, Industrial Park, ZhangSanZhai Town, ChangYuan City, XinXiang

City, HeNan province.

Name and address of Registered Manufacturer:

Yu An (Henan) Holding Limited Company

No. 168, Industrial Park, ZhangSanZhai Town, ChangYuan City, XinXiang

City, HeNan province.

Product name: Disposable Medical Protective Clothing

Product types: 160, 165, 170, 175, 180, 185

Product trademark: yūnn 字安控股

This certificate confirms that the product meets the requirements of the following standard

 The conformity was demonstrated based on following standard(s) listed by European Commission as harmonized with Directive 93/42/EEC.

EC Declaration of Conformity (Annex Class I according Rule 1 of Annex EN 14126:2003

VII of Directive 93/42/EEC) IX of Directive 93/42/EEC

The certification has been carried out In accordance with Individual rules and conditions agreed with the applicant. Evaluation has been carried out in accordance with:

Test report(s) No.: BG2003TR9184S06

Test conducted by: Shenzhen Bta Produce Testing Co.,Ltd.

Certificate issue date: 03.04.2020 Expiration date: 02.04.2025

NOTE:

- This certificate refers to the above mentioned produce and its conformity in regards of above mentioned standard(s) was proven on test sample
- This certificate does not Imply meeting all essential requirements, assessment of the series-production or any other restricted Notifled Bodies conformity assessment procedure appropriate for the product
- This certificate holder shall use this certificate in connection to declaration of conformity and technical data relevant for the product the certificate was Issued



ICR Polska Co.Ltd.

Plac Przymierza 6 03-944 Warsaw www.icrpolska.com e-mail: <u>icrpolska@icrqa.com</u>



Director: Rafal Kalinowski

Warsaw, 03. 04. 2020.





● CE 认证





● 宇安控股营业执照



● 宇安控股医疗器械生产许可证



● 宇安控股海关备案

海关进出口货物收发货人备案回执

企业名称	宇安 (河南) 控股有限公司
统一社会信用代码	91410728MA481B680B
海关备案日期	2020-04-16
海关编码	41079619AK
检验检疫备案号	4161100093
有效期	长期



自然人、法人或者非法人组织可通过"中国海关企业进出口信用信息公示平台"(http://credit.customs.gov.cn)或者"互联网+海关"(http://online.customs.gov.cn)查询海关公示的企业信息。

● 医疗器械生产产品登记表



▶ 一次性使用医用口罩 医疗器械注册证

中华人民共和国医疗器械注册证

注册证编号: 豫械注准 20202140676

注册证编	号: 豫械汪准 20202140070
注册人名称	宇安 (河南) 控股有限公司
注册人住所	河南省新乡市长垣市张三寨镇工业园 168 号
生产地址	河南省新乡市长垣市张三寨镇工业园 168 号
代理人名称	不适用
代理人住所	不适用的
产品名称	次性使用医用口罩
型号、规格	型号: 系带型、挂耳型; 规格: 大号、中号、小号。
结构及组成	产品由口罩体、鼻夹、口罩带组成,其中口罩体由无纺布及聚 丙烯熔喷布制成,鼻夹由可塑性材料制成,口罩带由弹性带或 无纺布制成。
适用范围	供临床各类人员在非有创操作过程中佩戴,覆盖住使用者的口鼻及下颌,为防止病原体微生物、颗粒物等的直接透过提供一定物理屏障。
附件	产品技术要求
其他内容	无
备注	22 日 16
+ 11 +0 2-1	The date of the first term of

审批部门: 河南省药品监督管理局

批准日期: 五〇二〇年三月九日 有效期至: 二〇二一年三月八日

● 一次性使用医用口罩 省级检验报告



No.202000680

河南省医疗器械检验所

检验报告



检验类别:应急检验

委 托 方: 宇安(河南) 控股有限公司

● 一次性使用医用口罩 合格证样张



产品名称 Product name	一次性使用医用口罩 Disposable Medical Face Mask
医疗器械生产许可证号 Medical Device Production License No.	豫药监械生产许20200216号/Henan Pharmaceutical Supervision Production License No.20200216
医疗器械注册证号 Medical Device Registration Certificate No.	豫械注准20202140676/Henan Pharmaceutical Supervision Standard 20202140676
执行标准 Executive Standard	YY/T 0969-2013 YY/T 0969-2013
型号 Model	挂耳型/Hanging Ear Type
规格 Specifications	175mm×95mm ±5% 175mm×95mm ±5%
数 <u>量</u> Quantity	50只/50 Masks
生产日期 Production Date	※ 按 服 、
生产批号 Production batch NO.	会 图 114
检验日期 Inspection Date	ZZ A
保质期 Quality Guarantee Period	2年/Two Years
检验员代号 Inspector Code	后於专用章

宇安(河南)控股有限公司 地址:中国·河南省新乡市长垣市张三寨镇工业园168号

YU'AN (HENAN) HOLDING LIMITED COMPANY

No.168,Zhangsanzhai Town Industrial Park,Changyuan City, Henan Province,China

ADD. No. 168, industrial park, zhangsanzhai Town, Changyuan City, Xinxiang City, Henan Province

● 宇安控股医用口罩 EN14683 检测报告

Report No.: BG2003TR9184S06 - Page 1 of 17 -



TEST REPORT EN 14683:2019 Medical face masks - Requirements and test methods Report Number...... BG2003TR9184S06 Tested by (+ signature) Roi Wei Jack chen Compiled by (+ signature) Jack Chen APPROVED Approved by (+ signature) Levi Hou Date of issue...... Mar. 20, 2020 Total number of pages...... 17 pages Testing Laboratory Shenzhen Bory Technology Service Co., Ltd GuangZhou Branch Address: 25F, Building 31, Changfeng International Zone, Xintang, Zengcheng District, Guangzhou Address As above Applicant's name Yu An (Henan) Holding Limited Company No. 168, Industrial Park, ZhangSanZhai Town, ChangYuan City, Address: XinXiang City, HeNan Province Test specification: Standard EN 14683:2019 Test procedure.....: Non-standard test method...... N/A Test Report Form No. EN 14683: 2019 Test Report Form(s) Originator: BORY Master TRF...... Dated 2020-02 The test results presented in this report relate only to the object tested. This report shall not be reproduced except in full, without the written approval of the Laboratory. The authenticity of this Test Report and its contents can be verified by contacting the Laboratory, responsible for this Test Report. Test item description...... Disposable Medical Mask Trade Mark Manufacturer Yang Ling Ao Rui Medical Devices Co., Ltd C5, Free Trade Street, Yangling Demonstration Zone, ShanXi Province, China. Model/Type reference I(clinical), II(surgical)

This Test Report is issued by the Company subject to its General Conditions of Service printed overleaf. Attention is drawn to the limitations of liability, indemnification and jurisdictional policies defined therein. The results shown in this test report refer only to the sample(s) tested unless otherwise stated and the sample(s) are retained for 30 days only.

Report No.: BG2003TR9184S06 - Page 2 of 17 -



List of Attachments (including a total number of pages in each attachment):

-- Attachment 1: One pages for Photo documentation.

Summary of testing:

Tests performed (name of test and test clause):

-- EN 14683:2019

Medical face masks - Requirements and test methods

Testing location:

25F, Building 31, Changfeng International Zone, Xintang, Zengcheng District, Guangzhou

Copy of marking plate

The artwork below may be only a draft. The use of certification marks on a product must be authorized by the respective NCBs that own these marks.

(Additional requirements for markings. See 1.7 NOTE)



Disposable Medical Mask

I (clinical)

EN 14683: 2019



Type

- Page 3 of 17 -

Report No.: BG2003TR9184S06

Report No.: BG2003TR9184S06 - Page 15 of 17 -



	EN 14683:2019		
Clause	Requirement Test	Result-Remark	Verdict
	Test specimens are complete masks or shall be cut from masks.	Complete mask.	Р
	Each specimen shall be able to provide 5 different circular test areas of 2,5 cm in diameter.	Diameter: 2,5 cm	Р
	If one specimen cannot provide 5 test areas of 2,5 cm in diameter, the number of test areas retrieved should be representative for the entire mask.	Diameter: 2,5 cm	Р
	The number of specimens that shall be tested is minimum 5, but can be greater and shall be increased if necessary to allow for an AQL of 4 %.	5 specimens	Р
	All specimens tested shall be taken from areas representative from the mask to incorporate all/any variation in construction.		Р
	Each test specimen shall be conditioned at (21 ±5) C and (85 ±5) % relative humidity for the time required to bring them into equilibrium with atmosphere prior to testing.	(21 ±2) °C; (85 ±2) %	Р
C.4	Procedure		Р
C.4.1	The test specimen is placed across the 2,5 cm diameter orifice (total area 4,9 cm2) and clamped into place so as to minimise air leaks and that the tested area of the specimen will be in line and across the flow of air.	Area: 4,9 cm2	P
C.4.2	The pump is started and the flow of air adjusted to 8 l/min.		Р
C.4.3	The manometers M1 and M2 are read and recorded.		Р
C.4.4	The procedure described in steps C.4.1 through C.4.3 is carried out on 5 (or appropriate number of) different areas of the mask and the readings averaged.		Р
C.5	Calculation of differential pressure		Р
	For each test specimen calculate the differential pressure ΔP as.	$\Delta P = (Xm1 - Xm2)/4,9$	Р
C.6	Test report	a number and date of this European Standard; b lot number or batch	Р
		code of the masks tested; c flow rate during	
		testing; d differential pressure for each test specimen.	

Report No.: BG2003TR9184S06 - Page 16 of 17 -



	EN 14683:2019		
Clause	Requirement Test	Result-Remark	Verdict
Ciausc	requirement rest	result-remain	Verdice
ZA	Relationship between this European Standard and the E EU Directive 93/42/EEC on medical devices	ssential Requirements of	Р
	This European Standard has been prepared under a mandate given to CEN by the European Commission Union to provide a means of conforming to the essential requirements of New Approach EU Directive 93/42/EEC concerning medical devices.		Р
	Once this standard is cited in the Official Journal of the European Union under that Directive and has been implemented as a national standard in at least one Member State, compliance with the clauses of this standard given in Table ZA.1 confers, within the limits of the scope of this standard, a presumption of conformity with the corresponding Essential Requirements of that Directive and associated EFTA regulations.		Р
ZA.1	Correspondence between this European Standard and E concerning medical devices	U Directive 93/42/EEC	Р
	Clause/subclause of this European Standard; 5.1.1, 5.1.2, 5.2.1, 5.2.2, 5.2.3, 6	Corresponding Essential Requirement of Directive 93/42/EEC: 8.1	Р
	Clause/subclause of this European Standard; 5.2.2	Corresponding Essential Requirement of Directive 93/42/EEC: 9.2	Р
	Clause/subclause of this European Standard; 6.	Corresponding Essential Requirement of Directive 93/42/EEC: 13	Р
	WARNING 1 — Presumption of conformity stays valid only as long as a reference to this European		Р
	Standard is maintained in the list published in the Official Journal of the European Union. Users of this		
	standard should consult frequently the latest list published in the Official Journal of the European		
	Union.		
	WARNING 2 —Other requirements and other EU Directives may be applicable to the product(s) falling within the scope of this standard.		Р

Report No.: BG2003TR9184S06 - Page 17 of 17 -



Photo-documentation





====== End of test report======

● 一次性使用医用口罩符合性声明

EC Declaration of Conformity

Manufacturer:

whose single Authorized EU-Representative:

Yu an (Henan) Holding Limited Company No. 168, Industrial Park, zhangsanzhai Town, Changyuan City, Xinxiang City, Henan Province, China

Luxus Lebenswelt GmbH
Kochstr.1, 47877, Willich, Germany
DIMID: DE/0000047791
Lin Sun
Tel: 0049- 1715605732
E-mail: info.m@luxuslw.de

We, the manufacturer, herewith declare that the products

Disposable Medical Mask

meet the provisions of Directive 93/42/EEC which apply to them.

The medical device has been assigned to class I according to Annex IX of the Directive 93/42/EEC. It bears the mark

((

following the procedure relating to the EC Declaration of Conformity set out in Annex VII of Directive 93/42/EEC.

This Declaration of Conformity covers all medical devices as specified in the product list belonging to this declaration.

The above mentioned declaration of conformity is exclusively under the responsibility of

Yu an (Henan) Holding Limited Company
No. 168, Industrial Park, zhangsanzhai Town, Changyuan City, Xinxiang City, Henan
Province, China

Place, date

Legally binding signature, Function

Changynan City May 6, 7000

MO7281025160

EC Declaration of Conformity

Page 1/1

● 医用外科口罩 医疗器械注册证

中华人民共和国医疗器械注册证

注册证编号: 豫械注准 20202140677

499-7-7-4	4. 4d. hd. 177. /r = 45.05.1.100.1.
注册人名称	宇安 (河南) 控股有限公司
注册人住所	河南省新乡市长垣市张三寨镇工业园 168 号
生产地址	河南省新乡市长垣市张三寨镇工业园 168 号
代理人名称	不适用
代理人住所	不适用
产品名称	医用外科口罩
型号、规格	型号: 挂耳型、系带型 规格: 175mm×95mm、175mm×90mm、140mm×90mm
结构及组成	产品分三层,外层由纺粘非织造布构成,中间层由聚丙烯熔喷布构成,另配有鼻夹、口罩带。
适用范围	用于戴在手术室医务人员口鼻部位,以防止皮屑、呼吸道微生物传播到开放的手术创面,并阻止手术病人的体液向医务人员传播,起到双向生物防护的作用。
附件	产品技术要求
其他内容	无
备 注	
审批部门;	河南省药品监督管理局 批准日期 一〇一〇年二日二十日

甲批部门:河南省药品监督管理局

批准日期: 二〇二〇年三月三十日 有效期至:二〇二一年三月二十九日 (审批部(T盖章)

● 医用一次性防护服 医疗器械注册证

注册人名称	宇安 (河南) 控股有限公司
	河南省新乡市长垣市张三寨镇工业园 168 号
	河南省新乡市长垣市张三寨镇工业园 168 号
代理人名称	The National Control of the Control
代理人住所	不 義用 曾 25
产品名称	医 用一次性防护服
型号、规格	型号: 连身式、分身式 规格: 160cm、165cm、170cm、175cm、180cm、185cm
结构及组成	防护服由帽子、上衣、裤子组成。
适用范围	用于医疗机构医护人员穿的职业防护衣,阻止来自患者的病毒 随空气或者液体向医务人员传播。
附件	产品技术要求
其他内容	无心心容然
备注	The state of the s
审批部门:	河南省药品监督管理局 批准日期3 二〇二〇年以月十五日 有效规章(一〇二十年以月十四日 (市北部门董章)

● 宇安控股工厂厂区图





● 医用一次性防护服 符合性说明

EC Declaration of conformity

2016/425 Personal protective equipment (PPE)

We,

Yu an (Henan) Holding Limited Company

No. 168, Industrial Park, Zhangsanzhai Town, Changyuan City, Xinxiang City, Henan

Province · China

Product Name: Disposable Medical Protective Clothing Product Model: 160, 165, 170, 175, 180, 185

The product has been assessed by the application of the following standards PPF:

EN 14126:2003+AC:2004, EN ISO 13982 -1:2004+A1:2010

Issue place and date

Company stamp and Signature of authorized personnel

Ohangynan City

Jiang Grang Sheng

May 6, 2020

GM